

### **Qualification of PET Scanner**

All imaging for ACRIN 6671 must be performed on a qualified PET instrument. Qualification instructions and application information can be found on the ACRIN Web site on the protocol specific page for ACRIN 6671 under “Imaging Materials”.

### **Qualification of MR Imaging**

Each site will be required to submit a test data set, of the four mandatory series, according to the trial specific protocol. The protocol parameter chart can be found on the ACRIN Web site, on the protocol specific page for ACRIN 6671 under “Imaging Materials”. Anonymized DICOM datasets for qualification must be sent on CD-ROM to the address below.

### **Removal of Confidential Participant Information**

The header record on DICOM formatted image data, which often contains information identifying the participant by name, MUST be scrubbed before the image is transferred. This involves replacing the Participant Name tag with the ACRIN Institution ID or number, replacing Participant ID tag with the ACRIN case number and putting the study number into the Other Participant ID tag. This can be performed using a customized software program available at your site, or TRIAD -- a software program available from ACRIN. Contact [Triad-Support@phila.acr.org](mailto:Triad-Support@phila.acr.org) with any questions.

### **Imaging Transmittal Worksheet**

An Imaging Transmittal Worksheet (ITW) must accompany all trial images submitted to ACRIN. The ITW is available on the ACRIN Web site on the protocol specific page for ACRIN 6671 under “Imaging Materials”. Images cannot be credited without a completed ITW.

### **Mail Images to:**

#### **American College of Radiology**

Core Laboratory  
1818 Market Street 16th floor  
Philadelphia, PA 19103  
Attn: 6671